

Troubleshooting

If problem cannot be corrected, consult your dealer/Homecare equipment provider.

PROBLEM	PROBABLE CAUSE	REMEDY
A. No flow	<ol style="list-style-type: none">1. Cylinder valve closed2. Conserver in "OFF" position3. Cylinder empty4. Dial not at setting5. Conserver not sensing breath	<ol style="list-style-type: none">1. Turn on cylinder2. Set to prescribed setting3. Replace cylinder4. Set dial so indicator points to a setting5. Check position of cannula in nose Do not breathe through mouth
B. Leaking at Yoke	<ol style="list-style-type: none">1. Missing or defective Yoke Seal	<ol style="list-style-type: none">1. Replace Yoke Seal

Precision Medical, Inc.
300 Held Drive • Northampton, PA 18067
800-272-7285 • 610-262-6090 • fax: 610-262-6080

504345rev3 25M 9/07 Printed in USA

Setting up your ***EASYPULSE5*** Oxygen Conserving Regulator



PRECISION  **MEDICAL**®



Quick Start Guide only
Refer to EasyPulse5 User Manual for complete operating instructions, warnings and cautions.


EasyPulse5 Oxygen Conserving Regulator (Conserver)

PARTS



STEP 1

Inspect Conserver and Post Valve for damage or contaminants. Do NOT use if any are found. "Crack" cylinder valve by opening the valve slightly and closing.



STEP 2

Ensure Conserver is in the OFF position and Yoke Seal is in place. Align Index Pins with Index Pin Alignment Holes in Post Valve. Hand tighten T-Handle. Do NOT use tools.



STEP 3

When using a Precision Medical Carry Bag, slide the cylinder into the bag with the gauge facing the mesh. Attach the cannula to the Conserver's Cannula Connection.



STEP 4

Slowly open the cylinder valve. Turn Dial to the prescribed setting. Zip the bag closed.



STEP 5

Place cannula in nostrils and breathe through nose to begin Oxygen flow. If using a bag adjust the straps for maximum comfort.



STEP 6

When finished using: Turn the cylinder valve to the OFF position. Turn Conserver Dial to "2 LPM Continuous" position. Wait for Oxygen to stop flowing from Conserver. Remove Conserver from Post Valve. Store Conserver and Carry Bag in a safe place.