










Mask Type	Nasal Pillows			Nasal Masks			Full Face Masks		
Suggested Customer Type	Should only be used with fixed pressure <8cmH ₂ O, patients on pressure >8cmH ₂ O might experience nasal irritation or nose bleeds. Not recommended for patients on blood thinning medication.			Suitable mask type for a majority of patients. Lightweight and comfortable as this mask type does not cover the whole face and promotes an effective breathing pattern with CPAP or Bi-level.			Recommended for patients with mouth leaks, blocked nasal passages, and deviated septum.		
Model	Swfit LT	Opus 360	OptiLife	Mirage Micro	Zest	ComfortGel	Mirage Quattro	Forma	ComfortFull
Diagram									
Range of Sizes	3 pillow sizes: small, medium, large	3 pillow sizes: small, medium, large	4 pillow sizes: petite, small, medium, large	5 sizes: small, medium, large, large-wide, extra-large	1 size; Easy-Clip, Premium Frosted silicone seal	4 sizes: petite, small, medium, large	4 sizes: extra-small, small, medium, large	4 sizes: small, medium, large, extra-large; Frosted silicone seal	3 sizes: small, medium, large
Cushions	dual-wall nasal pillow	built-in Soft-Lip cushion	soft, pillow cushion (removable)	dual-wall, with cushion clip	Foam FlexiFit cushion	customized ComfortGel™	Mirage™ dual-wall cushion	Foam FlexiFit cushion	customized ComfortGel™ (removable)
Elbow/Tube/Swivel	inter-changeable tube position	360° ball-and-socket elbow, inter-changeable tube position	inter-changeable tube position	360° elbow, quick release swivel	front moving swivel	360° swivel	simplified elbow with anti-asphyxia valve	front moving swivel	tubing quickly disconnected from the mask, flexible tubing
Forehead Support	N/A	N/A	N/A	"MicroFit dial" for a choice of 24 positions, providing a comfortable, effective fit	"Fit and Forget" permanent locking design	adjustable forehead support arm provides additional flexibility to meet the fitting needs of individual patients	"MicroFit dial" for a choice of 24 positions, providing a comfortable, effective fit	T-Piece	StabilitySelector forehead adjustment
Headgear	soft buckled headgear allows for easy adjustments, headgear stability arms, rotating barrel	three-dimensional support structures with hidden strap for quick release	two-coloured fabric headgear, two sets of headgear slots, chin support band supports seal	hypoallergenic Breath O-Prene™, 4 velcro fasteners, with 2 clips	Stretchgear™ headgear with Glider™ mechanism	soft flexible headgear straps improve comfort for all sleepers	set and forget headgear, 4 velcro fasteners with 2 clips	crown strap design with Glider™ mechanism, easy-release strap	premium headgear with EZ Peel™ tabs, updated ball-and-socket headgear attachments
Washing Care	hand wash recommended in warm soapy water	hand wash recommended in warm soapy water	hand wash the mask in warm water with a mild soap or dishwashing detergent	hand wash recommended in warm soapy water	hand wash recommended in warm soapy water	hand wash in warm water with mild dishwashing detergent, rinse thoroughly, and air dry	hand wash recommended in warm soapy water	hand wash recommended in warm soapy water	hand wash the mask in warm water with a mild dishwashing detergent